


The EKG Patch




Accurate. Fast. Consistent. Simple.

- ▶ Pre-positioned electrodes standardize application and minimize prep time.
- ▶ Built-in circuitry prevents lead wire reversal and placement errors.
- ▶ Single cable connection to existing EKG machines and portable devices.
- ▶ Latex-free, disposable materials promote, infection-free testing.
- ▶ Single application can remain in place for multiple readings.
- ▶ Immediate computerized interpretation of the EKG is available to unlimited resources.

The EKG Patch supports
Mason-Likar
electrode placement.

EKG Report



EKGs can be viewed
and stored on any IOS,
android, or PC device.

The Glasgow Interpretative
Algorithm provides a
detailed analysis of the
recorded EKG.

Your EKG Report can be
automatically uploaded
for Physician review and
interpretation.

Diagnostic, 12-Lead Resting EKGs in minutes... every time.

**TridentCare is partnering with MG Medical Products to offer the EKG Patch. The EKG Patch is an MG Medical Products product, and in order to receive the EKG Patch TridentCare customers will contract directly with MG Medical Products.*

The EKG Patch



All-in-One Design

Built-in lead wires attached to pre-positioned electrodes help ensure accurate lead placement on the chest wall and eliminates lead wire inversion errors. The EKG Patch comes in multiple sizes for optimal placement.

Flexible

One cable connects The EKG Patch to EKG devices in almost any care setting especially where time-efficient patient triage is critical.

Mitigate Unnecessary Hospital Admissions

The EKG Patch allows for on-demand EKG results in minutes, at any time and any place.

Minimal Technical Experience Required

FDA allows for application of the EKG Patch by anyone with little technical training. Single wire cable connects the EKG Patch to EKG devices, reducing potential for infection.



FDA-Cleared
CE-Marked
ISO 10993

www.tridentcare.com